

Assessing Technology Best Suited for Regenerative Medicine Biomanufacturing Biomonitoring


Aubrey Berger, Biomedical Engineering
Mentor: Dr. Vincent Pizziconi, Associate Professor
School of Biological and Health Systems Engineering

Introduction

- Regenerative medicine (RM) is an emerging and game-changing field of medicine that uses living cells as therapeutics in replaces or treating diseased tissue.
- Production of cells for the purpose of RM use is a complicated process that requires careful control of cell environments' pH, temperature, nutrient supply, waste disposal, sterility, and more.
- Special focus on stem cells, which are able to turn into many different kinds of cells.

Biomanufacturing:

- Bioreactors are the technological heart of biomanufacturing, providing controlled environment required for cell survival, proliferation and directed development of personalized RM treatments.



- Monitoring is a critical aspect in ensuring cells remain healthy and ideally characterized.
- FDA encourages use of process analysis technology (PAT)
 - A quality assurance system integrated into manufacturing process
 - Goal of ensuring final product quality through real time monitoring of critical quality factors, such as cell health, nutrient and waste levels, and other possible secretions¹.

Biomonitoring Techniques Considered

Dyes:


- Permanent and cytotoxic², even popular GFP

Imaging:

- Lack of ways to characterize cell other than count and vitality.
- Limited to 2D imaging


Scaffold Impedance Sensing:

- Incorporates cellular scaffold as sensor
 - Ohm's law: Voltage = Current * Impedance,
 - Capacitance Impedance = Capacitance * Rate of change of Voltage
- Impedance and capacitance change under different current frequencies and different cellular conditions
- Create an "Impedance or Capacitance Spectrum" to deduce cellular conditions⁴.



Impedance magnitude and phase over current frequency with and without cells and compared with different growth times⁴.

- Alternatively: track changes in impedance/capacitance over time to monitor changes in cell characteristics, such as the morphological transition from fibroblastic- to cuobidal-shaped cells as bones cells differentiate⁵.




Impedance magnitude over culture time with cultures of differentiating osteocytes in blues and control cells in red⁵.


- Relies on database of cellular characteristics correlated to their impedance/ capacitance spectrum, which is still under development.

Infrared Spectrophotometry:

- Beer-Lambert law: $A = \epsilon lc$
 - A = absorbance of light (at specific wavelength); ϵ = absorptivity, a property specific to material; l = pathlength of light beam through sample; c = concentration of substance
 - Absorption is proportional to concentration




- Measure absorption of light passing through cell culture with spectrophotometry



- Light beams, comprised of photon particles of various energy and momentum, can be differentiated by their wavelength
- Spectrophotometers, analytical instruments that use diffraction gratings to separate light beams into specific wavelengths, can be used to interrogate samples
- Measuring light absorbance at different wavelengths, allows multiple analytes to be measured at once in real time.
- Light near visible and infrared (IR) wavelengths especially promising for monitoring analytes of interest in biomanufacturing RM products:
 - Partial Least Square modeling of 1000nm-2500nm light proved accurate in predicting cell viability, glucose, lactate, and glutamate concentration *in vitro*⁶.
- At-line, real time monitoring of multiple parameters at once is possible, but requires advanced signal processing techniques.

Proposed IR Spectrophotometry Test System





Key Biomanufacturing Biomarkers

Glucose:

- During somatic reprogramming (turning a normal cell into a stem cell) there is a burst of glycolysis and glucose uptake⁷
- Naïve stem (not ready to differentiate) cells also consume more glucose than primed stem (ready to differentiate) cells.

Lactate:

- In stem cell cultures, lactate production differentiates between pluripotent stem stages: naive cells produce more than primed cells⁹



Design of Biomonitoring Test System

Next steps:

- Acquire spectrometer test system components
- Assemble and validate custom spectrometer system
- Test accuracy, reliability, rapidness, and feasibility of measuring glucose and lactate concentrations over IR range initially performed on mock solutions followed by testing in a bioreactor-like setting.
- Further development would be to test these same parameters across different scales of production and monitoring development of stem cells between naive and primed states.

References

[1] Food and Administration, D., Guidance for Industry PAT—A Framework for Innovative Pharmaceutical Development, Manufacturing, and Quality Assurance, DHHS, Rockville, MD 2004.
 [2] Dwarshuis, N. J., Parratt, K., Santiago-Miranda, A., & Roy, K. (2017). Cells as advanced therapeutics: State-of-the-art, challenges, and opportunities in large scale biomanufacturing of high-quality cells for adoptive immunotherapies. *Advanced Drug Delivery Reviews*, 114, 222-239. doi:10.1016/j.addr.2017.06.005
 [4] Tonello, Sarah et al. "Impedance-Based Monitoring of Mesenchymal Stromal Cell Three-Dimensional Proliferation Using Aerosol Jet Printed Sensors: A Tissue Engineering Application." *Materials* 13.10 (2020): 2231. Crossref. Web. <http://dx.doi.org/10.3390/ma13102231>.
 [5] Hildebrandt, B. (2010). Detection of the osteogenic differentiation of mesenchymal stem cells in 2D and 3D cultures by electrochemical impedance spectroscopy. *Journal of Biotechnology*, 148(1), 83-90. https://doi.org/10.1016/j.jbiotec.2010.01.007
 [6] Li, Mengyao et al. "Parallel Comparison of in Situ Raman and NIR Spectroscopies to Simultaneously Measure Multiple Variables Toward Real-Time Monitoring of CHO Cell Bioreactor Cultures." *Biochemical engineering journal* 137 (2018): 205-213. Web.
 [7] Ishida, T., Nakao, S., Ueyama, T. et al. Metabolic remodeling during somatic cell reprogramming to induced pluripotent stem cells: involvement of hypoxia-inducible factor 1. *Inflamm Regen* 40, 8 (2020). https://doi.org/10.1186/s41232-020-00117-8
 [8] Spiegelman, Clifford H et al. "Theoretical Justification of Wavelength Selection in PLS Calibration: Development of a New Algorithm." *Analytical chemistry* (Washington) 70.1 (1998): 35-44. Web.
 [9] Tsogtbaatar, Enkhkuul et al. "Energy Metabolism Regulates Stem Cell Pluripotency." *Frontiers in cell and developmental biology* 8 (2020): 87-87. Web.
 [10] Mamouei, Mohammad et al. "Comparison of Wavelength Selection Methods for in-Vitro Estimation of Lactate: a New Unconstrained, Genetic Algorithm-Based Wavelength Selection." *Scientific reports* 10.1 (2020): 16905-16905. Web.

Image Credits

Sunglasses by Andy Mc from the Noun Project
 sun by tulpahn from the Noun Project
 lightbulb by Maxim Kulikov from the Noun Project
 Flask by Alex Muravev from the Noun Project
 camcorder by Graphic Engineer from the Noun Project
 EM spectrum from Philip Roman
 Lens by Hoang Loi from the Noun Project
 Light Bulb by Jacob from the Noun Project
 BioReactor by Emily Earll from the Noun Project
 Bioreactor by f Dr. Yu-Shik Hwang
 Cells by Adnen Kadri from the Noun Project